



Ear Pressure Relief Device

- Ideal for children and adults
- Opens the Eustachian Tube with air
- Delivers a constant airflow into the nasal cavity
- For patients with problems associated with:
 - Flying
 - Allergies and Colds
 - Scuba Diving
 - Swimming



Professional Version
EP-3000



Home Version
for patients (Rx only)
EP-2100



Tips for Professional Version
EP-1000

The EarPopper® Device is a safe, simple device used to open the Eustachian Tube and equalize pressure in the middle ear. This device is the non-invasive, non-drug related way to treat Otitis Media with Effusion, Eustachian Tube Dysfunction and other issues related to middle ear pressure. Ideal for children and adults.

The EarPopper delivers a constant flow of air into the nasal cavity and then the patient swallows. During the moment of swallowing, air is diverted up the Eustachian Tube, ventilating the middle ear. The EarPopper instantly relieves negative ear pressure and allows accumulated fluids to drain.

Treatment can be performed in the office with the Professional Version (EP-3000) and prescribed for continued treatment at home with the Home Version (EP-2100).

Product Number	Description	Units / Box
EP-3000	EarPopper, Professional Version <i>includes 10 disposable tips and 1-9 volt battery</i>	1
EP-2100	EarPopper, Home Version (<i>for patients, Rx in US only</i>) <i>includes 4-AAA batteries</i>	1
EP-1000	Disposable tips	50

For treatment guidelines or patient information please visit
www.earpopper.com

To order call **1.888.229.2875**



Summit Medical, Inc. | St. Paul, MN USA | local p: 651.789.3939 f: 651.789.3979 | toll free p: 1.888.229.2875 f: 1.888.229.1941

www.summitmedicalusa.com

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.
©2013 Summit Medical, Inc. EarPopper is a registered trademark of Summit Medical, Inc. All rights reserved. EP-0010/E JAN13