

Shikani Middle Meatal Antrostomy Stent



- Maintains patency of the middle meatal antrostomy
- Removes in the office without anesthesia
- Reduces or eliminates the need for post-operative packing
- Soft, medical-grade silicone is easily compressed and manipulated to facilitate insertion and maximize patient comfort
- Eliminates adhesions caused by lateralization of the middle turbinate



Developed and refined in use at Johns Hopkins Hospital by Dr. Alan Shikani, this stent is designed to be placed following functional endoscopic sinus surgery. Maintain patency of antrostomy and minimize adhesions between the middle turbinate and lateral nasal wall with this proven stent.^{1,2}

Product Number	Description	Unit/Box <i>Sterile, single use</i>	Dimension Details (in millimeters)																													
SP-80010	Small	6		<table border="1"> <thead> <tr> <th>Measurement</th> <th>Small</th> <th>Medium</th> <th>Large</th> </tr> </thead> <tbody> <tr> <td>A</td> <td>27.9</td> <td>35.7</td> <td>39.0</td> </tr> <tr> <td>B</td> <td>5.0</td> <td>6.0</td> <td>8.1</td> </tr> <tr> <td>C</td> <td>21.9</td> <td>29.5</td> <td>29.5</td> </tr> <tr> <td>D</td> <td>6.2</td> <td>7.8</td> <td>9.3</td> </tr> <tr> <td>E</td> <td>4.6</td> <td>4.6</td> <td>4.6</td> </tr> <tr> <td>F</td> <td>23.4</td> <td>33.4</td> <td>33.4</td> </tr> </tbody> </table>	Measurement	Small	Medium	Large	A	27.9	35.7	39.0	B	5.0	6.0	8.1	C	21.9	29.5	29.5	D	6.2	7.8	9.3	E	4.6	4.6	4.6	F	23.4	33.4	33.4
Measurement	Small	Medium		Large																												
A	27.9	35.7		39.0																												
B	5.0	6.0	8.1																													
C	21.9	29.5	29.5																													
D	6.2	7.8	9.3																													
E	4.6	4.6	4.6																													
F	23.4	33.4	33.4																													
SP-80020	Medium	6																														
SP-80030	Large	6																														

Literature and Support Materials:

¹ Shikani, A., *A New Middle Meatal Antrostomy Stent for Functional Endoscopic Sinus Surgery*. Laryngoscope, 104:638:641, 1994

Doyle, D., *Middle Meatal Spacer and Septal Splint, Operative Techniques in Otolaryngology-Head and Neck Surgery*, Vol. 8, No 2 (JUN), 1997

² Shikani, A., *Middle Antrostomy Stenting Following Pediatric Endoscopic Sinus Surgery*. Am J Rhinol, 10:225-228, 1996 – CD ROM.

Shikani, A., *Clinical Use of the Shikani Middle Meatal Antrostomy Stent in the Pediatric and Adult Patients* – CD ROM.

Call **1.888.229.2875** to order today!



Summit Medical, Inc. | St. Paul, MN USA | local p: 651.789.3939 f: 651.789.3979 | toll free p: 1.888.229.2875 f: 1.888.229.1941

www.summitmedicalusa.com

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.
©2011 Summit Medical, Inc. All rights reserved. ENT-0014/B APR11