



# Foreign Body Removal

Designed to safely and easily remove objects from ear and nasal cavities with suction.




**CPT 69200\***  
(Ear) Removal foreign body from external auditory canal; without general anesthesia

**CPT 30300\***  
(Nasal) Removal foreign body, intranasal; office type procedure

## Vacutract™

Foreign Body Removal Device



- Hooks up to standard suction tubing
- For office, ER or any area with access to a suction system

Product #	Description	Units / Box
VA-1100	Vacutract (with tip) <i>Non-sterile, single use</i>	5

## VacuTIP™

Foreign Body Suction Tip



- Suction cup conforms to almost any shape
- Disposable tips can be used on most standard 7 French Baron Suction

Product #	Description	Units
VA-1120	VacuTIP <i>Non-sterile, single use</i>	20 tips / case

### RELATED PRODUCTS

Product Number	Description	Units / Box
DV-7305D-D	DeVilbiss® Suction Pump, AC only	1
DV-7305P-D	DeVilbiss® Suction Pump, AC and Battery	1
SI-62-22-20	Baron Suction, 7 French	1



\*The information contained in this document is provided for convenience only and represents no statement, promise, or guarantee by Summit Medical, Inc. concerning coverage or levels of reimbursement. It is always the provider's responsibility to determine coding, coverage and claim information for services provided. For CY 2012 Program Fee Schedules through Centers for Medicare & Medicaid Services (CMS) visit [www.cms.gov](http://www.cms.gov). Please note: fee schedule corrections and charges occur periodically. CPT copyright 2012 American Medical Association. All rights reserved. CPT is a registered trademark of the American Medical Association. Applicable FARS/DFARS restrictions apply to government use.

To order call **1.888.229.2875**



Summit Medical, Inc. | St. Paul, MN USA | local p: 651.789.3939 f: 651.789.3979 | toll free p: 1.888.229.2875 f: 1.888.229.1941

[www.summitmedicalusa.com](http://www.summitmedicalusa.com)

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. ©2013 Summit Medical, Inc. Vacutract and VacuTIP are trademarks of Summit Medical, Inc. All rights reserved. ENT-0018/E DEC13